Parcus 3.5mm Hex Driver with Handle

Important Product Information

Produktinformationsblatt

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Parcus 3.5mm Hex Driver with Handle

1. Indications:

The Parcus 3.5mm Hex Driver with Handle is designed to insert and remove Parcus Titanium Interference Screws. The driver features a 1.5mm cannulation and is compatible with Parcus 1.5mm Guide Wires.

2. Warnings:

- A. This product is intended for use by or on the order of a physician.
- B. This product should not be used to implant any device other than Parcus Titanium Interference Screws.
- C. Do not use this product if it is bent or there is visible damage to the shaft or drive section.

3. Material:

This device is comprised of a Stainless Steel Alloy shaft and an Aluminum handle. The materials used in the manufacture of this device that are intended to be placed inside the body are radio-opaque and can, therefore, be detected with conventional X-Ray or fluoroscopy.

4. Inspection:

- A. Inspect the device for damage at all stages of handling.
- B. If damage is detected, consult the manufacturer for guidance.

5. Cleaning:

- A. Immediate rinsing and cleaning after use with an enzymatic detergent will effectively remove and prevent drying of adherent blood, tissue, etc.
- B. Scrub device with a soft brush, paying special attention to the cannulation and other areas where debris might accumulate. Always avoid any harsh materials that can scratch or mar the surface of the device.
- C. Rinse the device thoroughly with water following the cleaning process.

6. Sterilization:

This device is not provided sterile. It should be cleaned and sterilized per facility standards using validated protocols.

7. Directions for use:

- A. Place a Parcus Titanium Interference Screw onto the drive portion of the driver until it is fully engaged.
- B. Place the screw and driver over an appropriately placed Parcus 1.5mm Guide Wire to guide the screw into the insertion site.
- C. Rotate the driver clockwise until the screw has reached the desired depth. A counterclockwise rotation of the driver will remove the screw.
- D. Remove the driver from the screw, followed by removal of the guide wire.

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