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Parcus 35 Knotless PEEK CF Suture Anchor Driver

1. Indications:

Parcus 35 Knotless PEEK CF Suture Anchor Drivers are provided sterile as complements to 35 Knotless PEEK CF Suture Anchors. They are used to assist in the insertion of 35 Knotless PEEK CF Suture Anchors in bone.

2. Warnings:

- A. This product is intended for use by or on the order of a physician.
- B. This product is designed for single use only and should not be re-sterilized. Drivers that are reused can fatigue and break.
- C. Never use a bent driver.
- D. This product is not intended for use as an implant.
- E. Loading the driver in an off axis manner may result in breakage of the implant or bending of the driver.

3. Material:

This device is comprised of a Stainless Steel shaft and with a Plastic handle. The shaft of this device is manufactured from a Stainless Steel Alloy and is

radio-opaque. The shaft can, therefore, be detected with conventional X-Ray or fluoroscopy.

4. Sterilization:

This device is provided in sterile packaging. The contents are sterilized by EO gas. This product must never be re-sterilized.

5. Packaging and Labeling:

- A. Do not use this product if the packaging or labeling has been damaged, shows signs of exposure to moisture or extreme temperature or has been altered in any way.
- B. Please contact Parcus Medical Customer Service to report any package damage or alterations.

6. Storage:

Products must be stored in the original unopened package in a dry place and must not be used beyond the expiration date indicated on the package.

7. Directions for use:

- A. Aseptically open the package containing the 35 Knotless PEEK CF Suture Anchor and the package containing the 35 Knotless PEEK CF Suture Anchor Driver.
- B. Assemble the suture anchor to the driver by inserting the threaded tip of the driver into the socket of the anchor and rotating the driver in a

clockwise direction until the driver is firmly attached to the anchor.

Warning: Off-axis or insufficient seating of the anchor to the driver may result in breakage of the implant.

- C. Pass sutures and insert implant into bone using the preferred user technique as described in implant DFU.
- D. Once implant is fully inserted, remove driver by rotating the driver in a counterclockwise direction until the driver is free from the anchor.

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