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4.5mm V-LoX™ PEEK CF Hybrid Suture Anchor

1. Indications

The Parcus 4.5mm V-LoX PEEK CF Hybrid Suture Anchor is indicated for attachment of soft tissue to bone. This product is intended for the following indications:

Shoulder Rotator Cuff Repair,
Acromioclavicular Separation Repair,
Bankart Lesion Repair, Biceps
Tenodesis, Capsular Shift or
Capsulolabral Reconstruction, Deltoid
Repair, SLAP Lesion Repair.

Knee Medial Collateral Ligament Repair,
Lateral Collateral Ligament Repair,
Posterior Oblique Ligament Repair,
Extra Capsular Reconstruction,
Iliotibial Band Tenodesis, Patellar
Ligament and Tendon Avulsion
Repair.

Foot/Ankle Lateral Stabilization, Medial
Stabilization, Midfoot Reconstruction,
Achilles Tendon Repair, Hallux
Valgus Reconstruction, Metatarsal
Ligament Repair.

Elbow Tennis Elbow Repair, Biceps Tendon Reattachment.

Hand/Wrist Scapholunate Ligament Reconstruction, Ulnar or Radial Collateral Ligament Reconstruction, TFCC.

2. Contraindications

- A. Any active infection.
- B. Blood supply limitations or other systemic conditions that may retard healing.
- C. Foreign body sensitivity, if suspected should be identified and precautions observed.
- D. Insufficient quality or quantity of bone. Suture anchor performance is directly related to the quality of bone into which the anchor is placed.
- E. Patient's inability or unwillingness to follow the surgeon's prescribed post-operative regimen.
- F. Any situation that would compromise the ability of the user to follow the directions for use or using the device for an indication other than those listed.

3. Adverse Effects

- A. Infection, both deep and superficial.
- B. Allergies and other reactions to device materials.
- C. Risks due to anesthesia.

4. Warnings

- A. This product should only be used by or on the order of a physician.
- B. The fixation provided by this device should be protected until healing is complete. Failure to follow the postoperative regimen prescribed by the surgeon could result in the failure of the device and compromised results.
- C. Anchor size selection should be made with care taking into consideration the quality of the bone in which the anchor is to be placed.
- D. Any decision to remove the device should take into consideration the potential risk of a second surgical procedure. Adequate postoperative management should be followed after implant removal.
- E. Polyethylene polyblend sutures, as used in this device, elicit a minimal inflammatory reaction in tissues, followed by gradual encapsulation of the suture by fibrous connective tissue. Polyethylene polyblend suture is not absorbed, nor is any significant change in tensile strength retention known to occur in vivo.
- F. The patient should be advised of the use and limitations of this device.
- G. Pre-operative planning and evaluation, surgical approaches and technique, and familiarity of

the implant, including its instrumentation and limitations are necessary components in achieving a good surgical result.

- H. This device must never be reused.
- I. This device must never be re-sterilized.
- J. Appropriate instrumentation should be used to implant this device.

5. Packaging and Labeling

- A. Do not use this product if the packaging or labeling has been damaged, shows signs of exposure to moisture or extreme temperature or has been altered in any way.
- B. Please contact Parcus Medical Customer Service to report any package damage or alterations.

6. Material Specifications

The 4.5mm PEEK CF Push-In Hybrid Suture Anchor is manufactured from PEEK CF (Polyetheretherketone carbon fiber reinforced) and is supplied with two Polyethylene polyblend sutures. In knotless applications, as described below, these sutures are removed after the anchor has been driven into position. A disposable stainless steel driver with molded plastic handle is included.

7. Sterilization

The 4.5mm V-LoX PEEK CF Hybrid Suture Anchor is supplied in sterile packaging. The contents are sterilized by EO gas. This product must never be re-sterilized.

8. Storage

This product must be stored in the original unopened package in a dry place and must not be used beyond the expiration date indicated on the package.

9. Directions for Use (Knotless Application)

- A. Identify bone of sufficient quantity and quality into which the anchor is to be placed.
- B. The Parcus 4.5mm V-LoX PEEK CF Hybrid Suture Anchor is inserted into a socket that is created by a Drill Bit or Awl/Punch. The standard Drill Bit or Awl/Punch is 3.2mm. A 4.0mm Drill Bit is used in extremely dense bone. Arthroscopic cannulae are typically used in conjunction with these devices. Always try to approach the targeted site as close to perpendicular as possible.
- C. Place sutures through the targeted tissue. These sutures will be secured when the anchor is driven into the prepared socket. The 4.5mm V-LoX PEEK CF Hybrid Suture Anchor will accommodate up to (6) strands of #2 USP suture. In arthroscopic applications, draw the

- sutures up through the same cannula through which the implant will be used.
- D. If using the Drill Bit, insert the Drill Bit into a drill chuck or pin driver to a depth that covers the laser etched line on the proximal end of the Drill Bit, as applicable. Advance the Drill Bit into the bone approximately 16mm or until the circumferential laser mark on the distal tip of Drill Bit is flush with the surface of the bone, as applicable. Remove the Drill Bit. If using the Awl/Punch, advance its tip to the distal circumferential laser mark. As an option, the Drill Bit or Awl/Punch may be left in position until ready to insert the anchor.
 - E. Aseptically open the package containing the 4.5mm V-LoX PEEK CF Hybrid Suture Anchor.
 - F. Pass the sutures placed through the tissue described in Step C through the eyelet in the distal tip of the suture anchor. If the Drill Bit or Awl was left in position as described in Step D, remove the Drill Bit or Awl before proceeding. The proper orientation of the suture and implant is to have the suture pass perpendicular to the plane of the tissue through which it is passed with the hole of the eyelet aligned with the path of the suture.

- G. Place the distal end of the suture anchor into the previously created hole and advance by gently tapping until the circumferential laser mark on the distal tip of the driver is flush with the bone. **Failure to insert the driver to the depth indicated by the laser mark may result in the anchor being left proud.**
- H. Remove the two #2 sutures provided with the anchor and discard or use for other means of fixation. Cut the tails of the remaining sutures.
- I. Remove all remaining instrumentation.

10. Directions for Use (Knotted Application)

- A. Identify bone of sufficient quantity and quality into which the anchor is to be placed.
- B. The Parcus 4.5mm V-LoX PEEK CF Hybrid Suture Anchor is inserted into a socket that is created by a 3.2mm drill bit or Awl/Punch. Arthroscopic cannulae are typically used in conjunction with these devices. Always try to approach the targeted site as close to perpendicular as possible.
- C. If using the Drill Bit, insert the Drill Bit into a drill chuck or pin driver to a depth that covers the laser etched line on the proximal end of the Drill Bit, as applicable. Advance the Drill Bit into the bone approximately 16mm or until the circumferential laser mark on the distal tip of Drill Bit is flush with the surface of the bone, as

applicable. Remove the Drill Bit. If using the Awl/Punch, advance its tip to the distal circumferential laser mark. As an option, the Drill Bit or Awl/Punch may be left in position until ready to insert the anchor.

- D. Aseptically open the package containing the 4.5mm V-LoX PEEK CF Hybrid Suture Anchor.
- E. If the Drill Bit was left in position as described in Step C, remove the Drill Bit before proceeding. Place the distal end of the suture anchor into the previously created hole and advance by gently tapping. Once the anchor is started into the hole, you may continue to drive implant in or you may rotate driver in a clockwise direction and thread implant into bone until the circumferential laser mark on the distal tip of the driver is flush with the bone.

Warning: Insertion of the anchor off axis from the hole may result in implant failure. Failure to insert the driver to the depth indicated by the laser mark may result in the anchor being left proud.

- F. Pass the two #2 sutures provided with the implant through tissue and tie down. Cut the tails of the sutures just above the knots.
- G. Remove all remaining instrumentation.

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