## Parcus Femoral Offset Guide

### **Important Product Information**

Produktinformationsblatt

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## **Parcus Femoral Offset Guide**

### **1. Indications:**

The Parcus Femoral Offset Guides are used to assist in placing a 2.4mm guide pin for creating the femoral tunnel or socket in ACL reconstructions.

### **2. Warnings:**

- A. This product is intended for use by or on the order of a physician.
- B. No component of this instrument is intended to be left in the body.
- C. Never use this device if it has been bent or damaged in any way.
- D. Extending the knee with this instrument in position may result in damage to the instrument.
- E. Take care to remove all debris from the cannulation of the instrument after use.

### **3. Material:**

This product is manufactured from stainless steel and aluminum alloys. The materials used in the manufacture of this device that are intended to be placed inside the body are radio-opaque and can, therefore, be detected with conventional X-Ray or fluoroscopy.

**4. Sterilization and Cleaning:**

This device is not provided sterile. It should be cleaned and sterilized per facility standards using validated protocols.

**5. Directions for use:**

- A. Select the proper offset. The Femoral Offset Guides are sized in various millimeter offsets that place a guide pin anterior to the tang of the guide. For example, a 7mm aimer will leave a 2mm wall of bone if a 10mm reamer is used.
- B. Insert the aimer into the knee through a tibial tunnel or arthroscopic portal. It may be necessary to flex the knee to place the guide in the desired position.
- C. Place the tang of the Femoral Offset Guide at the “over the top” position on the femur.
- D. Advance a 2.4mm guide pin through the cannulation in the guide into or through the bone.
- E. Remove the Femoral Offset Guide and use a cannulated drill or reamer over the guide pin to create the tunnel or socket in the targeted bone.

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