Parcus Graft Fixation System (GFS)

Important Product Information

Produktinformationsblatt

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Parcus Graft Fixation System (GFS)

1. Indications:

The Parcus GFS is indicated for use in the fixation of ligaments and tendons in patients requiring ligament or tendon repair.

2. Contraindications:

- A. Any active infection.
- B. Blood supply limitations or other systemic conditions that may retard healing.
- C. Foreign body sensitivity, if suspected, should be identified and precautions observed.
- D. Insufficient quality or quantity of bone.
- E. Patient's inability or unwillingness to follow surgeon's prescribed post-operative regimen.
- F. Any situation that would compromise the ability of the user to follow the directions for use or using the device for an indication other than those listed.

3. Adverse Effects:

- A. Infection, both deep and superficial.
- B. Allergies and other reactions to device materials.
- C. Risks due to anesthesia.

4. Warnings:

- A. This product is intended for use by or on the order of a physician.
- B. The fixation provided by this device should be protected until healing is complete. Failure to follow the postoperative regimen prescribed by the surgeon could result in the failure of the device and compromised results.
- C. Size selection of the implant should be made with care taking into consideration the quality of the bone

on which the implant will rest and the desired length of graft that will reside in the tunnel.

- D. Any decision to remove the device should take into consideration the potential risk of a second surgical procedure. An adequate postoperative management plan should be implemented after implant removal.
- E. Pre-operative planning and evaluation, surgical approaches and technique, and familiarity of the implant, including its instrumentation and limitations are necessary components in achieving a good surgical result.
- F. This device must never be reused.
- G. This device must never be re-sterilized.
- H. Appropriate instrumentation should be used to implant this device.

5. Packaging and Labeling:

- A. Do not use this product if the packaging or labeling has been damaged, shows signs of exposure to moisture or extreme temperature or has been altered in any way.
- B. Please contact Parcus Medical Customer Service to report any package damage or alterations.

6. Material Specifications:

The Parcus GFS is supplied with high-strength, braided, polyethylene polyblend sutures. The anchor material is Ti-6Al-4V ELI (ASTM F136).

7. Sterilization:

This Parcus GFS is supplied sterile.

8. Storage:

Products must be stored in the original unopened package in a dry place and must not be used beyond the expiration date indicated on the package.

9. Directions for Use (Soft Tissue Graft)

- A. Create a bone tunnel, sized to accommodate the graft, using the preferred technique.


Note: The Parcus GFS is supplied in sizes to be used with specific tunnel diameters. See package label for details.

- B. Measure the overall length of the tunnel with the Parcus GFS Depth Gauge.
- C. Subtract the amount of graft to be placed into the tunnel from the overall tunnel length.
- D. The remainder will be the length of the loop required to position the desired amount of graft in the tunnel. If “T” equals the tunnel length and “G” equals the desired length of graft to be inserted into the tunnel and “L” equals the length of Loop required to leave the desired graft in the tunnel then $T-G=L$. For example:

Tunnel length =	45mm
Graft in Tunnel =	25mm
Loop length=	20mm

- E. Aseptically open the appropriately sized Parcus GFS. Remove the foil wrapper from the suture loops and draw one end of the prepared soft tissue graft, through the space between the titanium implant and the loops of sutures. Parcus recommends whip-stitching the free ends of the graft to aid in tensioning the graft (see step “K”).
- F. Position the loops of suture at the midpoint of the length of the graft.
- G. Orient the GFS/Graft construct in such a way that the axis of the titanium implant is parallel to the graft. Place the wedge shaped end of the GFS at the zero point on the GFS Depth Gauge and make a mark on

the graft with a surgical marker that corresponds to the length of the femoral tunnel described as dimension “T” above. This will indicate the point at which the wedge shaped end of the GFS has passed the lateral femoral cortex. See Drawing “A” below.

- H. When using the Parcus GFS for ACL or PCL reconstruction, select a suture passing device that will allow the user to pass the leading suture of the GFS through the tibial tunnel and exit out the femoral tunnel (e.g. Parcus Drill Tip Guide Pins with Eyelet).
- I. Once the leading suture can be retrieved, grasp it with a hemostat or similar instrument and pull the graft into the femoral tunnel until the mark, described in “G” passes by the tunnel opening.
- J. Using the sutures attached to the free ends of the graft, pull the graft in a caudad direction. This movement will allow the wedge shaped end of the GFS to catch on the femoral cortex and orient the GFS over the tunnel opening. Pull vigorously on the ends of the graft to determine if the GFS is properly deployed. If fixation is not achieved remove the graft and repeat step “I”.
- K. Once secure fixation has been verified, remove the leading suture.
- L. Tension the graft by pulling the sutures.
- M. Fix the distal end of the graft in the tibial tunnel according to user preference.

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