## Parcus Threaded Dilator

### Important Product Information

Produktinformationsblatt

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## **Parcus Threaded Dilator**

### **1. Indications:**

The Parcus Threaded Dilator is designed to create a starter hole to ease the insertion of Interference Screws in Bone-Tendon/Ligament-Bone ACL Reconstructions. The Threaded Dilator is cannulated and is compatible with 1.5mm or less Guide Wires.

### **2. Warnings:**

- A. This product is intended for use by or on the order of a physician.
- B. Do not use this product if it is bent or there is visible damage to the shaft or drive section.
- C. This device is only to be used when sufficient bone quality and implant/graft size exists such that fixation is not compromised.
- D. This device is not intended to be used as an implant.

### **3. Material:**

This device is comprised of a Stainless Steel Alloy shaft and an Aluminum handle. The materials used in the manufacture of this device that are intended to be placed inside the body are radio-opaque and can, therefore, be detected with conventional X-Ray or fluoroscopy.

**4. Inspection:**

- A. Inspect the device for damage at all stages of handling.
- B. If damage is detected, consult the manufacturer for guidance.

**5. Cleaning:**

- A. Immediate rinsing and cleaning after use with an enzymatic detergent will effectively remove and prevent drying of adherent blood, tissue, etc.
- B. Scrub device with a soft brush, paying special attention to areas where debris might accumulate. Always avoid any harsh materials that can scratch or mar the surface of the device.
- C. Rinse the device thoroughly with water following the cleaning process.

**6. Sterilization:**

This device is not provided sterile. It should be cleaned and sterilized per facility standards using validated protocols.

**7. Directions for use:**

- A. Place the Threaded Dilator over an appropriately placed 1.5mm Guide Wire to guide the dilator into the insertion site.
- B. Rotate the Threaded Dilator clockwise until the dilator has reached the desired depth.

- C. Rotate the Threaded Dilator counterclockwise to remove the dilator.
- D. Place the desired Interference Screw over the Guide Wire and insert per applicable Directions for Use.

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