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Parcus Drill Bits for Suture Anchor Placement

1. Indications:

Parcus Drill Bits are provided in various lengths and diameters as complements to suture anchors and/or Drill Guides. They are used to create sockets and tunnels in bone.

2. Warnings:

- A. This product is intended for use by or on the order of a physician.
- B. This product is designed for single use only and should not be re-sterilized. Drill Bits that are reused can fatigue and break.
- C. Never use a bent Drill Bit.
- D. This product is not intended for use as an implant.
- E. Loading the Drill Bit in an off axis manner may result in breakage or bending.

3. Material:

This device is manufactured from a Stainless Steel Alloy. The materials used in the manufacture of this device that are intended to be placed inside the body are radio-opaque and can, therefore, be detected with conventional X-Ray or fluoroscopy.

4. Inspection:

- A. Inspect the device for damage at all stages of handling.
- B. If damage is detected, consult the manufacturer for guidance.

5. Sterilization:

This device is not provided sterile. It should be cleaned and sterilized per facility standards using validated protocols.

6. Packaging and Labeling:

- A. Do not use this product if the packaging or labeling has been damaged, shows signs of exposure to moisture or extreme temperature or has been altered in any way.
- B. Please contact Parcus Medical Customer Service to report any package damage or alterations.

7. Storage:

Products must be stored in the original unopened package in a dry place and must not be used beyond the expiration date indicated on the package.

8. Directions for use:

- A. Place the Drill Bit into the chuck of a purpose designed drill or wire driver. Drill Bits may have laser etched lines to help determine depth of penetration. Attach Drill Bits so that laser lines are visible as required.
- B. If a Drill Guide is utilized, advance the Drill Bit through the guide and into the bone to the required depth.
- C. Remove the Drill Bit from the guide.
- D. Remove the Drill Guide, as necessary.

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