## Parcus Awl/Punch

### Important Product Information

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## Parcus Awl/Punch

### 1. Indications:

The Parcus Awl/Punch is designed to gauge patient bone quality and/or to create sockets of specific size in bone.

### 2. Warnings:

- A. This product is intended for use by or on the order of a physician.
- B. Do not use this product if it is bent or there is visible damage to the shaft or drive section.
- C. This device is not intended to be used as an implant.

### 3. Material:

This device is manufactured from a Stainless Steel Alloy and may include an aluminum alloy handle. The materials used in the manufacture of this device are radio-opaque and can, therefore, be detected with conventional X-Ray or fluoroscopy.

### 4. Inspection:

- A. Inspect the device for damage at all stages of handling.
- B. If damage is detected, consult the manufacturer for guidance.

### 5. Cleaning:

- A. Immediate rinsing and cleaning after use with an enzymatic detergent will effectively remove and prevent drying of adherent blood, tissue, etc.
- B. Scrub device with a soft brush, paying special attention to areas where debris might accumulate. Always avoid any harsh materials that can scratch or mar the surface of the device.
- C. Rinse the device thoroughly with water following the cleaning process.

### 6. Sterilization:

- A. This device is not provided sterile. It should be cleaned and sterilized using validated protocols before each use.
- B. The following table provides recommended minimum sterilization parameters that have been validated by Parcus Medical to provide a  $10^{-6}$  sterility assurance level (SAL):

Cycle Type	Minimum Temperature	Minimum Exposure Time (Wrapped)	Minimum Dry Time
Gravity	132°C / 270°F	15 Minutes	30 Minutes <sup>1</sup>
Pre-Vacuum	132°C / 270°F	4 Minutes	

<sup>1</sup> Drying times vary according to load size and should be increased for larger loads.

### 7. Instructions for use:

- A. Place the tip of the Awl/Punch on the surface of the bone at the point at which the bone quality is to be gauged or the socket is to be created.
- B. Use a suitable mallet to insert the Awl/Punch into the bone to the desired depth or until the circumferential laser line is flush with the bone.
- C. Remove the Awl/Punch.

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