**Parcus Single Use, Unsterile Drill Bits
for Suture Anchor Placement**

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www.ParcusMedical.com
(941)755-7965

0086
EC REP
MDSS GmbH
Schiffgraben 41
30175 Hanover, Germany

Parcus Drill Bits for Suture Anchor Placement

1. Indications:

Parcus Drill Bits are provided in various lengths and diameters as complements to suture anchors and/or Drill Guides. They are used to create sockets and tunnels in bone.

2. Warnings:

- A. This product is intended for use by or on the order of a physician.
- B. This product is designed for single use only and should not be re-sterilized. Drill Bits that are reused can fatigue and break.
- C. Never use a bent Drill Bit.
- D. This product is not intended for use as an implant.
- E. Loading the Drill Bit in an off axis manner may result in breakage or bending.

3. Material:

This device is manufactured from a Stainless Steel Alloy. The materials used in the manufacture of this device that are intended to be placed inside the body are radio-opaque and can, therefore, be detected with conventional X-Ray or fluoroscopy.

4. Inspection:

- A. Inspect the device for damage at all stages of handling.
- B. If damage is detected, consult the manufacturer for guidance.

5. Cleaning:

- A. Immediate rinsing and cleaning after use with an enzymatic detergent will effectively remove and prevent drying of adherent blood, tissue, etc.
- B. Scrub device with a soft brush, paying special attention to areas where debris might accumulate. Always avoid any harsh materials that can scratch or mar the surface of the device.
- C. Rinse the device thoroughly with water following the cleaning process.

6. Sterilization:

This device is not provided sterile. The following table provides recommended minimum sterilization parameters that have been validated by Parcus Medical to provide a 10^{-6} sterility assurance level (SAL):

Cycle Type	Minimum Temperature	Minimum Exposure Time (Wrapped)	Minimum Dry Time
Gravity	132°C / 270°F	15 minutes	30 minutes ¹
Pre-vacuum	132°C / 270°F	4 minutes	

¹ Drying times vary according to load size and should be increased for larger loads.

7. Packaging and Labeling:

- A. Do not use this product if the packaging or labeling has been damaged, shows signs of exposure to moisture or extreme temperature or has been altered in any way.
- B. Please contact Parcus Medical Customer Service to report any package damage or alterations.

8. Storage:

Products must be stored in the original unopened package in a dry place and must not be used beyond the expiration date indicated on the package.

9. Directions for use:

- A. Place the Drill Bit into the chuck of a purpose designed drill or wire driver. Drill Bits may have laser etched lines to help determine depth of penetration. Attach Drill Bits so that laser lines are visible as required.
- B. If a Drill Guide is utilized, advance the Drill Bit through the guide and into the bone to the required depth.
- C. Remove the Drill Bit from the guide.
- D. Remove the Drill Guide, as necessary.

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