



PEEK CF Interference Screws are recommended for the attachment of soft tissue to bone.

Features & Benefits

Made of carbon fiber-reinforced PEEK-OPTIMA

- Non-absorbable, radiolucent, and MR safe
- Modulus of elasticity closely matching cortical bone
- Naturally lubricious, minimizes “graft clocking”

Fully-threaded design

- Optimal fixation for both soft tissue and BTB grafts

Rounded edge

- Protects soft tissue grafts

Driver is cannulated, tapered & square

- Provides a positive driver-to-screw interface
- Distributes insertion torque and eliminates “driver capture”

PEEK CF Interference Screws

Part #	Diameter (mm)	Length (mm)	Description	Material	QTY/Box
10021	7	20	Interference screw	PEEK CF	1
10022	7	25	Interference screw	PEEK CF	1
10023	7	30	Interference screw	PEEK CF	1
10025	8	20	Interference screw	PEEK CF	1
10026	8	25	Interference screw	PEEK CF	1
10027	8	30	Interference screw	PEEK CF	1
10028	8	35	Interference screw	PEEK CF	1
10029	9	20	Interference screw	PEEK CF	1
10030	9	25	Interference screw	PEEK CF	1
10031	9	30	Interference screw	PEEK CF	1
10032	9	35	Interference screw	PEEK CF	1
10034	10	25	Interference screw	PEEK CF	1
10035	10	30	Interference screw	PEEK CF	1
10036	10	35	Interference screw	PEEK CF	1
10166	11	25	Interference screw	PEEK CF	1
10167	11	30	Interference screw	PEEK CF	1
10168	11	35	Interference screw	PEEK CF	1
10038	12	25	Interference screw	PEEK CF	1
10039	12	30	Interference screw	PEEK CF	1
10040	12	35	Interference screw	PEEK CF	1

Instrumentation for PEEK CF Interference Screws

Part #	Diameter (mm)	Joint	Description	Sterile/Unsterile	Single-use/Reusable
10838	4.3	Knee	Tunnel notcher	Unsterile	Reusable
10275	-	Knee	PEEK CF interference screw driver shaft, w/1.5mm cannulation	Unsterile	Reusable
10042	-	Knee	PEEK CF interference screw driver shaft, hudson, w/1.5mm cannulation	Unsterile	Reusable
10372	6.4	Knee	Interference screw threaded dilator, w/1.5mm cannulation	Unsterile	Reusable



Knee