



GFS II is a suspensory fixation device recommended for use in the fixation of ligament or tendon repairs. The larger button may eliminate the stepped tunnel technique required by smaller buttons. It is an excellent choice for revision surgeries.

## Features & Benefits

### Loop made of Ultra High Molecular Weight Polyethylene (UHMWPE)

- High strength
- Non-absorbable

### Button made of 6Al-4V ELI titanium alloy

- Provides the strength and stiffness required for successful ligament reconstruction

### Pre-loaded with color-coded pulling suture and flipping suture

- Convenience in the operating room

GFS II Large					
Part #	Width of Button (mm)	Length of Button (mm)	Description	Length of Loop (mm)	QTY/Box
10702	4.4	16.5	GFS, II large	12	1
10703	4.4	16.5	GFS, II large	15	1
10704	4.4	16.5	GFS, II large	20	1
10705	4.4	16.5	GFS, II large	25	1

Instrumentation for GFS II Large					
Part #	Diameter (mm)	Joint	Description	Sterile/Unsterile	Single-use/Reusable
10057	-	Knee	GFS Depth Gauge	Unsterile	Reusable



**Knee**