



Twist & Triple Twist PEEK Screw-In Suture Anchors are recommended for use in both large and small-joint repairs.

## Features & Benefits

### Fully-threaded design

- Optimal cortical fixation
- Reduce the risk of "pull-back"

### Made of PEEK-OPTIMA Natural

- Non-absorbable, radiolucent, MR safe

#### Twist PEEK Screw-In Suture Anchors

Part #	Diameter (mm)	Length (mm)	Description	Material	QTY/Box
10511T	4.5	12.6	Twist, screw-in, suture anchor, w/2, #2 sutures (blu & wht/grn)	PEEK	1
10517T	5.5	12.6	Twist, screw-in, suture anchor, w/2, #2 sutures (blu & wht/grn)	PEEK	1
10518T	6.5	12.6	Twist, screw-in, suture anchor, w/2, #2 sutures (blu & wht/grn)	PEEK	1
10521	4.5	12.6	Twist, screw-in, suture anchor, w/2, #2 sutures (wht & wht/blu), w/needles (MO-6)	PEEK	1
10522	5.5	12.6	Twist, screw-in, suture anchor, w/2, #2 sutures (wht & wht/blu), w/needles (MO-6)	PEEK	1
10523	6.5	12.6	Twist, screw-in, suture anchor, w/2, #2 sutures (wht & wht/blu), w/needles (MO-6)	PEEK	1

#### Triple Twist PEEK Screw-In Suture Anchors

Part #	Diameter (mm)	Length (mm)	Description	Material	QTY/Box
10524T	5.5	12.6	Triple Twist, screw-in, suture anchor, w/3, #2 sutures (blu, wht/grn & wht/blk)	PEEK	1

#### Twist PEEK Screw-In Suture Anchors w/Two Suture Tapes

Part #	Diameter (mm)	Length (mm)	Description	Material	QTY/Box
11003	4.5	12.6	Twist, screw-in, suture anchor, w/2, 1.6mm suture tapes (wht/blu & wht/blk)	PEEK	1
11283	4.5	12.6	Twist, screw-in, suture anchor, w/2, 1.6mm suture tapes (blu, wht/blu)	PEEK	1
11285	4.5	12.6	Twist, screw-in, suture anchor, w/2, 1.6mm suture tapes (blk, wht/blk)	PEEK	1
11004	5.5	12.6	Twist, screw-in, suture anchor, w/2, 1.6mm suture tapes (wht/blu & wht/blk)	PEEK	1
11284	5.5	12.6	Twist, screw-in, suture anchor, w/2, 1.6mm suture tapes (blu, wht/blu)	PEEK	1
11286	5.5	12.6	Twist, screw-in, suture anchor, w/2, 1.6mm suture tapes (blk, wht/blk)	PEEK	1
11005	6.5	12.6	Twist, screw-in, suture anchor, w/2, 1.6mm suture tapes (wht/blu & wht/blk)	PEEK	1

#### Twist SST PEEK Screw-In Suture Anchors w/One Suture Tape

Part #	Diameter (mm)	Length (mm)	Description	Material	QTY/Box
11349	4.5	12.6	Twist SST, screw-in, suture anchor, w/1, 2mm suture tape (blu)	PEEK	1
11350	4.5	12.6	Twist SST, screw-in, suture anchor, w/1, 2mm suture tape (blk)	PEEK	1
11351	5.5	12.6	Twist SST, screw-in, suture anchor, w/1, 2mm suture tape (blu)	PEEK	1
11352	5.5	12.6	Twist SST, screw-in, suture anchor, w/1, 2mm suture tape (blk)	PEEK	1

#### Instrumentation for Twist PEEK Screw-In Suture Anchors

Part #	Diameter (mm)	Joint	Description	Sterile/Unsterile	Single-use/Reusable
10365	3.2	Shoulder	Awl	Unsterile	Reusable
10338	4.1	Shoulder	Awl	Unsterile	Reusable
10609	4.5	Shoulder	Punch/tap	Unsterile	Reusable
10277	5.5	Shoulder	5.5mm or 6.5mm punch/tap	Unsterile	Reusable
10222	5.5	Shoulder	5.5mm or 6.5mm punch/tap, hudson	Unsterile	Reusable
10505U	3.2	Small Joint	Drill bit, w/positive stop	Unsterile	Reusable
10540U	3.2	Shoulder	Drill bit, w/positive stop	Unsterile	Reusable
10541U	4	Shoulder	Drill bit	Unsterile	Reusable
10823	-	Small Joint	8-point drill guide, w/beaver tail handle	Unsterile	Reusable
11000U	4	Small Joint	Drill bit	Unsterile	Reusable
11241	4	Small Joint	Drill guide	Unsterile	Reusable



**Shoulder**